Labormöbel  
Lab furniture

## LGT-QE-SERIE

## LGT-QE-SERIES

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**Bestellnummer:** LGT-QE-MH

## Contents

<b>1</b>	<b>Important safety instructions</b>	<b>13</b>
<b>2</b>	<b>Important notice to consumers</b>	<b>13</b>
<b>3</b>	<b>Operating conditions and conventional usage</b>	<b>13</b>
<b>4</b>	<b>Models</b>	<b>14</b>
<b>5</b>	<b>Technical data</b>	<b>14</b>
	Current consumption: .....	14
	Internal power outlet(s): .....	14
	Disconnecting device: .....	15
	Sound proofing: .....	15
<b>6</b>	<b>Scope of supply</b>	<b>15</b>
<b>7</b>	<b>Installation</b>	<b>15</b>
	Positioning .....	15
	Height compensation and locking of the castors .....	15
	Connection of the External Display Module APPS-EDU (optional).....	16
	Installation of the vacuum pump .....	16
<b>8</b>	<b>Function test</b>	<b>16</b>
<b>9</b>	<b>Maintenance</b>	<b>16</b>
<b>10</b>	<b>APPS</b>	<b>17</b>
<b>11</b>	<b>Accessories (not included)</b>	<b>19</b>
	External Display Unit for APPS-Systems .....	19
	Interface cable for PC-connection .....	19
	Monitor .....	19

## 1 Important safety instructions

- Read these instructions.
- Keep these instructions.
- Heed all warnings.
- Follow all instructions.
- Do not use this apparatus near water.
- Mains powered apparatus shall not be exposed to dripping or splashing.
- Do not block any ventilation openings. Install in accordance with the manufacturer's instructions.
- Do not defeat the safety purpose of the polarized or grounding-type plug. A polarized plug has two blades with one wider than the other. A grounding type plug has two blades and a third grounding prong. The wide blade or the third prong are provided for your safety. If the provided plug does not fit into your outlet, consult an electrician for replacement of the obsolete outlet. (for USA and Canada)
- Protect the power cord from being walked on or pinched particularly at plugs, convenience receptacles, and the point where they exit from the apparatus.
- **WARNING:** This apparatus shall be connected to a mains socket outlet with a protective earthing connection only.
- Refer all servicing to qualified service personnel. Servicing is required when the apparatus has been damaged in any way, such as power-supply cord or plug is damaged, liquid has been spilled or objects have fallen into the apparatus, the apparatus has been exposed to rain or moisture, does not operate normally, or has been dropped.
- Use only with accessories specified by the manufacturer, or sold with the apparatus.
- Use only power cables rated at least for the maximum input current of the apparatus.
- If the apparatus is used in a manner not specified by the manufacturer the protection provided by the equipment may be impaired.

## 2 Important notice to consumers

- This product has been manufactured according to strict specifications and voltage requirements that are applicable in the country in which it is intended that this product should be used. If you have purchased this product via the internet, through mail order, and/or via a telephone sale, you must verify that this product is intended to be used in the country in which you reside.
- **WARNING:** Use of this product in any country other than that for which it is intended could be dangerous and could invalidate the manufacturer's or distributor's warranty. Please also retain your receipt as proof of purchase otherwise your product may be disqualified from the manufacturer's or distributor's warranty.

## 3 Operating conditions and conventional usage

- The integrated noise reduction cabinet is exclusively designed as a noise reduction device for vacuum pumps.
- The vacuum pump may not have tangible dangerous parts.
- The bench must be installed on an even and horizontal surface.
- The incoming and exhausted airflow must be able to stream unhindered. The frontal air intakes and the rear air discharge openings must be distant from other objects at least 10 cm.
- The fans must be able to turn free and may not be blocked.
- The bench may only be used in clean and dry areas.
- The heat emission of the installed pump(s) may not exceed 1.5kW (models with 4 fans) respectively 2.5kW (models with 8 fans).
- The bench may not be operated above 30°C / 86°F ambient temperature and above 80% humidity.
- The inner temperature of the noise reduction cabinet may not exceed 65 °C / 149 °F
- During transport and storage temperature and humidity must be within the following range:  
Transport: Temperature: -20 to 60°C (-4 to 140°F), humidity: 5 to 90%  
Storage: Temperature: 10 to 40°C (50 to 104°F), humidity: 10 to 80%

- If condensation water is formed after unpacking the bench, an acclimatisation period of 2 hours must be adhered prior to installation.
- Easily inflammable and combustible liquids are not allowed inside the noise reduction cabinet.
- If a liquid should penetrate the electrical part of the bench it is to be set immediately out of operation and must be checked by a specialist.
- The valid accident guarding regulations according to each country must be considered.
- A repair of the bench may be accomplished only by a specialist.

## 4 Models

### Format of the model-numbering:

**Lab-Bench:** LGTQE-<W>x<D>x<H>

<W>: Width in cm

<D>: Depth in cm

<H>: Height in cm

**Noise Reduction Cabinet:** LGTQE-SS<P>-<F><C>U<N>

<P>: Mounting position

L: Left

R: Right

<F>: Number of built-in fans

4: Four fans for a maximum heat emission of 1.5kW

8: Eight built-in fans for a maximum heat emission of 2.5kW

<C>: Electrical connectors

1: IEC C13/C14 plug/receptacle. Max. current: 10A

2: IEC C19/C20 plug/receptacle. Max. current: 16A

<N>: Number of connections

1: Model for one pump -> One power inlet, one power outlet

2: Model for two pumps -> Two power inlets, two power outlets

## 5 Technical data

### Current consumption:

Model	$I_{max}$	$I_{oc}$
LGTQE-SSX-X1UX	10A	300mA
LGTQE-SSX-X2UX	16A	300mA

$I_{max}$  is the maximum current per power inlet that can be drawn with a load connected to the inner power outlet(s).  $I_{oc}$  is the maximum current drawn by the electronics without a load connected to the inner power outlet(s).

### Internal power outlet(s):

The internal power outlets are intended for the connection of electrical loads inside the noise reduction cabinet. The length of the power cord may not exceed 3m (9 ft). The current load may not exceed 10A/16A according to the model.

The voltage of the power outlet is the same as the line voltage the noise reduction cabinet is connected to.

### Disconnecting device:

The power cable(s) serves as the disconnecting device. To completely disconnect the equipment from the electric supply network disconnect all power cables. Please assure that the power cables are accessible at all times.

### Sound proofing:

Environmental friendly recycling product. Oil and water repellent. Inflammability fulfills FMVSS 302 and DIN 75200.

## 6 Scope of supply

- 1 x Lab Bench LGTQE-Series
- 1 x Power cord IEC10A respectively IEC16A (models with noise reduction cabinet)
- 1 x Extension cord IEC C13/C14 10A respectively IEC C19/C20 16A (models with noise reduction cabinet)

## 7 Installation

### Positioning

Position the table at the required place. If the table has to be lifted, do not grab it at the tabletop but use the frame. As the massive frame is quite heavy, the tabletop could be damaged. Benches with integrated noise reduction cabinet have to be placed in a way that the frontal air intakes and the rear air outtakes are not blocked. However the table can be placed directly to a wall as the tabletop overlaps and so the minimum distance to the wall is guaranteed.

### Height compensation and locking of the castors


Fig. 3: Adjustment of the castors

① Grub screw

② Plate

③ Break

If the floor is not perfectly even and for that reason not all castors have contact to the floor, they can be prolonged up to 20mm. To do so, please loosen the grub screw on the appropriate castor by approximately one turn and turn the plate clockwise to extend the table-leg. One complete turn extends the leg by 1.5mm. After the adjustment tighten the grub screw again. To ensure a safe stand, please lock the break of all castors.

## Connection of the External Display Module APPS-EDU (optional)

If the External Display Module APPS-EDU has been ordered as well, it can be connected to the table now. To do so, please connect it to the appropriate socket on the back of the noise reduction cabinet and place it on the tabletop. As soon as the table is connected to a power outlet all the data collected like exhaust temperature, speed of the fans, current consumption of the fans etc. will be shown on the display. Also the behaviour of the APPS-module can be adjusted here.

## Installation of the vacuum pump

Open the door of the noise reduction cabinet by turning the quick fasteners counterclockwise. Now completely pull out the oil tray and place the vacuum pump centric on it. The oil sight glass should face the front and the motor to the rear of the cabinet. Feed the vacuum- and exhaust hoses through the feedthroughs in the tabletop and/or the back panel. While doing so, please make sure that the hoses are long enough to connect the pump. The oil tray can now be slid back in. Avoid squeezing and crippling of any cables and hoses while doing so. Now connect the cables and hoses accordingly to the manual of your instrument and connect the pump to the inner power outlet. Please keep in mind that as soon as the table is connected to a power outlet, also the inner socket is powered, so please make sure to connect the table to the power no sooner as you want to power the pump as well. If your pump has a power plug that doesn't fit the power outlet, an adapter cable can be order from Sonation. Alternatively the power cable can be fed through one of the feedthroughs and connected to a power outlet directly. If installed that way, the noise reduction cabinet has to be connected to power outlet seperately to ensure the cooling of the installed pump.

## 8 Function test

After installation of the noise reduction box, the instrument can be started up again. Please check if the LED on the back panel is flashing every five seconds and the fans are running. To get an acoustical confirmation, please briefly push the button on the back panel. After 5 seconds at the latest a confirmation beep should sound. Depending on the chosen OpMode the box beeps one, two or three times. If another beepcode is audible, please refer to the APPS section in this manual. If the LED is not flashing and the fans are not running, please check the connection to the mains.

After 4 to 6 hours of operation the pump has reached its operating temperature. After this time please assure yourself again of the correct function of the box by checking that the LED flashes every 5 seconds and the fans are running. If the external display module APPS-EDU is used, the correct function of the electronics can be looked up on the display.

## 9 Maintenance

To clean the surfaces please use a soft and moist tissue.

### every day

A malfunction is indicated by an optical and acoustical alert. If an alert is showing up, please refer to the chapter APPS in this manual and follow the instructions accordingly to the alert code.

### every month

Check the air intakes, fans and protection grids for any kind of dirt. Clean as necessary. To clean the grids, the pump must be switched off and the noise reduction cabinet must be disconnected from the power source.

## 10 APPS

APPS stands for Active Pump Protection System and offers a comprehensive monitoring of all security-related functions of the pump and the noise reduction cabinet. All benches with noise reduction cabinet are equipped with this system by default.

Fig. 4: APPS-Module

- ① Interface
- ② Piezo-speaker
- ③ Push-button with LED
- ④ Reset

The Interface is used to connect to a PC or to the optional external display module.

**For the connection to a PC the optional interface cable APPS-ICBL is needed. Connection to ethernet devices or to a telephone may cause damage to the connected devices.**

### The following data and functions are monitored during operation:

- Exhaust temperature
- Function and speed of the fans
- Power supply and function of the noise reduction box
- Current consumption of the connected load

There are three user selectable modes. The modes can be switched by pressing the push button for at least 10 seconds until a confirmation signal sounds. Depending on the chosen mode the box beeps one, two or three times.

When critical values are exceeded an optical and acoustical warning shows up. The following list shows the different warning codes and their meaning.

By pressing the push-button ③ an actual alarm can be muted. Depending on the severity level it is muted for one hour or one day. The LED keeps blinking representing the beepcode.

**Blink/Beepcode: -                      Meaning: APPS module defective**

While connected to the mains supply the LED has to blink every 5 seconds for one, two or three times.

If not, please press the reset button ④ through the small hole next to the push-button ③ by using a thin screwdriver or a bent up paper clip. If this does not work, please contact Sonation to discuss further steps.

**Blink/Beepcode: 1                      Meaning: OpMode 1: Normal operation**

The noise reduction cabinet monitors all relevant functions. This mode is the factory default.

By pressing the push-button ③ the function of the piezo speaker can be checked. After 5 seconds at the latest the speaker should beep one time, confirming OpMode 1.

**Blink/Beepcode: 2                    Meaning: OpMode 2: Age-related wear monitoring deactivated**

In OpMode 1 the age related wear of the fans is monitored. When the fan speed falls permanently below 80% of the original fan speed an alert shows up. In OpMode 2 this alert is suppressed. Nevertheless on a total failure of the fans a different alert still shows up.

**Blink/Beepcode: 3                    Meaning: OpMode 3: Age related wear monitoring + emergency shutdown deactivated**

In OpMode 1 and OpMode 2 the internal power outlet is switched off if the temperature of the exhaust air reaches a critical value (default: 60°C) or if too much current is drawn by the load connected to the internal power outlet (default: 12A at 10A models, 19,2A at 16A models)

In OpMode 3 this emergency shutdown is deactivated.

**This mode may only be used when the expected damage of shutting down the pump will be higher than the damage that will occur by an overheated pump. When using this mode the user has to assure the proper function of the pump at any time.**

**Blink/Beepcode: 2-1                Meaning: Power fail**

The external power supply has failed.

If the mains supply has been broken intentionally, the alert can be muted by pressing the push-button ③. As soon as the capacitors are discharged, the module will shut down automatically. With the acoustical alert switched off this will happen after around 5 hours.

If the mains supply has not been broken intentionally, please check the fuses on the machine or the power outlet the cabinet is connected to.

**Blink/Beepcode: 2-2                Meaning: Temperature warning (default: T>40°C)**

The exhaust temperature is higher than the programmed alert limit. Make sure the air intakes and outlets are not blocked, the pump does not produce too much waste heat and the ambient temperature is not too high.

**Blink/Beepcode: 2-3                Meaning: Overtemperature (default: T>60°C)**

The exhaust temperature is higher than the programmed critical limit. The internal power outlet is shut down (exception: OpMode 3). Make sure the air intakes and outlets are not blocked, the pump does not produce too much waste heat and the ambient temperature is not too high. Check the proper function of the fans and the pump.

**Blink/Beepcode: 3-1                Meaning: Life expectancy of the fans reached (n < 80% n0)**

One time a week the maximum fan speed is tested automatically. If the motor speed falls below 80% of the original motor speed an alert shows up for about half a minute (exception: OpMode 2 and 3). The back panel should be exchanged to assure the proper cooling of the insulated pump.

**Blink/Beepcode: 3-2                Meaning: Fan failure or fan speed too low (n < 700 rpm)**

If the fan speed of one or more fans falls below 700 rpm this alert shows up. Make sure the fans are not blocked.

**Blink/Beepcode: 3-3                Meaning: Complete fan failure**

This error occurs if all fans fail. Open the front door of the noise reduction cabinet to prevent overheating of the insulated pump. The back panel should be replaced immediately.

**Blink/Beepcode: 1-2-1              Meaning: Configuration error**

The number of power supplies, temperature sensors and fans are pre-programmed to the module. This alert shows an erroneous configuration. Please contact Sonation to discuss further steps.

**Blink/Beepcode: 1-2-2              Meaning: Discharged or defective capacitors**

After connecting the box to the mains supply the backup capacitors are getting charged. During the charging procedure (about 10 minutes) this alert can show up. If the alert is still there after this time, please contact Sonation to discuss further steps.

**Blink/Beepcode: 1-2-3      Meaning: Failure of a temperature sensor**  
To protect the pump all fans run on 100% speed. The back panel should be replaced immediately.

**Blink/Beepcode: 1-3-1      Meaning: +24V voltage supply erroneous**  
The internal power supply is erroneous. The back panel has to be replaced immediately.

**Blink/Beepcode: 1-3-2      Meaning: +5V voltage supply erroneous**  
The internal voltage regulator is erroneous. The back panel has to be replaced immediately.

**Blink/Beepcode: 2-1-2      Meaning: Overcurrent (default: 12A/19,2A)**  
The power consumption of the internal power outlet is monitored. If the current consumption is higher than the programmed value for a longer time (~30seconds) the power outlet is switched off to protect the electronics and the connected pump. As factory default this value is set to 12A at 10A models and to 19,2A at 16A models.

## 11 Accessories (not included)

### External Display Unit for APPS-Systems

The display can be connected to all LGT-QE lab benches with noise reduction cabinet. All data that is monitored by the APPS-System is clearly represented on the screen.

**Order number:** APPS-EDU

### Interface cable for PC-connection

All LGT-QE lab benches with noise reduction cabinet can be connected to a PC to present and process the collected data with Sonation or third party software.

**Order number:** APPS-ICBL

### Monitor

Mount your monitor ergonomically, space-saving and securely with this height-adjustable swivel arm.

**Order number:** LGT-QE-MH

## IMPROVE YOUR SCIENCE

SONATION GMBH  
ALTE SCHULSTRASSE 39  
DE-88400 BIBERACH

TEL.: +49 (0) 7351 30038-39  
E-MAIL: MAIL@SONATION.COM  
WWW.SONATION.COM